WHY THE REMNOVA COMPLETE BITE BLOCK IS VALUABLE FOR NEUROANESTHESIA

With MEP monitoring, it is imperative to have a bite block to prevent complications from the repeated biting caused by the neuro-stimulation.

The *Complete Bite Block* has many advantages over the gauze bite block. To name a few:

- the soft plastic device is designed to fit securely in patients' mouths, ensuring proper bite block placement even when flipping prone (it is also compatible with the prone pillow)
- third party tested to withstand four times the maximum human bite force
- helps stabilize the endotracheal tube leading to a more secure airway
- allows access to the airway for intubation, suctioning, oral airway placement, etc.
- fits under the face mask to allow for bag mask ventilation
- many more on the website remnova.com - including links to the related studies

